# Uri-Cath<sup>TM</sup>

## Closed, Pre-assembled Urinary Drainage System

# Because babies in the NICU *require* gentle urinary drainage



Accurate measurement of urinary output for neonates and pediatric patients

Soft silicone material minimizes tissue irritation and encrustation

- Pre-assembled, sterile system reduces risk of infection and reduces valuable nursing time
- Stopcock with sampling access site allows for non-contaminated midstream specimen sampling
- 150 ml burette for accurate urine output measurement

UTAH MEDICAL



Utah Medical Products, Inc manufactures
a broad range of neonatal critical care
devices including the Gesco™ family of
products

PRODUCTS INC.



The Uri-Cath<sup>™</sup> tray consists of a pre-connected silicone urinary drainage catheter, three-way stopcock for sampling purposes, and a 150 ml burette for accurate urine output measurement.

#### **Design Features**

- Biocompatible silicone tubing
- Encased radiopaque stripe
- Soft, rounded tip
- Stopcock with latex-free injection site
- 150 ml collection burette with drainage stopcock
- 3 distal openings
- Available in four sizes

### **Ordering Information**

#### **Urinary Drainage Trays**

Model #	Description	Unit of Measure	Image
4193507	3.5 French Uri-Cath™ set with silicone catheter	10/BOX	
4195007	5.0 French Uri-Cath™ set with silicone catheter	10/BOX	
4196507	6.5 French Uri-Cath™ set with silicone catheter	10/BOX	
4198007	8.0 French Uri-Cath™ set with silicone catheter	10/BOX	

### **Urinary Drainage Catheters**

Model #	Description	Unit of Measure	Image
4193505	3.5 French Uri-Cath™ catheter only	10/BOX	
4195005	5.0 French Uri-Cath™ catheter only	10/BOX	
4196505	6.5 French Uri-Cath™ catheter only	10/BOX	
4198005	8.0 French Uri-Cath™ catheter only	10/BOX	

#### **United States:**

Utah Medical Products, Inc. 7043 South 300 West Midvale, Utah 84047 800.533.4984 • 801.566.1200 www.utahmed.com Europe: ECREP
Utah Medical Products Ltd.
Athlone Business & Technology Park
Athlone, Co. Westmeath N37 XK74
Ireland
353.90.647.3932

